

Tetramer Lab Core Standard Protocols

Instructions for shipping peptides:

Peptides should be provided either in lyophilized (dry powder) or in a reconstituted DMSO solution at 20 mg per ml. We require a minimum of 1 mg of dry powder (please provide exact weight of peptide) or 50 uL of DMSO solution to process an order.

Peptide in DMSO needs to be shipped to BRI in dry ice. Dry peptide can be sent in a 0.5 Eppendorf tube via regular FedEx (no dry ice).

Addendum:

If you need to order the peptide, greater than 80% should be good. We don't have a preferred company, but you can order that through Mimotopes, and have the peptide shipped directly to us.